

For Research Use Only

GMPS Recombinant antibody, PBS Only

Catalog Number: 86130-1-PBS



Basic Information

Catalog Number:

86130-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9350

GenBank Accession Number:

BC012178

GeneID (NCBI):

8833

UNIPROT ID:

P49915

Full Name:

guanine monphosphate synthetase

Calculated MW:

693 aa, 77 kDa

Observed MW:

67 kDa

Purification Method:

Protein A purification

CloneNo.:

250704D4

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86130-1-RR (GMPS antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86130-1-PBS in a different storage buffer formulation.