

For Research Use Only

RAD23A Recombinant antibody

Catalog Number: 86132-1-RR



Basic Information

Catalog Number: 86132-1-RR	GenBank Accession Number: BC014026	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 5886	CloneNo.: 250716C7
Source: Rabbit	UNIPROT ID: P54725	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: RAD23 homolog A (S. cerevisiae)	
Immunogen Catalog Number: AG0716	Calculated MW: 40 kDa	
	Observed MW: 55 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Jurkat cells, MCF-7 cells, SKOV-3 cells, OVCAR-3 cells, C6 cells, mouse testis tissue, rat testis cells
Species Specificity: human, mouse, rat	

Background Information

RAD23A, also named as HHR23A, belongs to the RAD23 family. It is a multiubiquitin chain receptor involved in modulation of proteasomal degradation. RAD23A binds to 'Lys-48'-linked polyubiquitin chains in a length-dependent manner and with a lower affinity to 'Lys-63'-linked polyubiquitin chains. It is proposed to be capable to bind simultaneously to the 26S proteasome and to polyubiquitinated substrates and to deliver ubiquitinated proteins to the proteasome. It is involved in nucleotide excision repair and is thought to be functional equivalent for RAD23B in global genome nucleotide excision repair (GG-NER) by association with XPC. In vitro, the XPC:RAD23A dimer has NER activity. Can stabilize XPC. It is also involved in vpr-dependent replication of HIV-1 in non-proliferating cells and primary macrophages and is required for the association of HIV-1 vpr with the host proteasome.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various cell lysates were subjected to SDS PAGE followed by western blot with 86132-1-RR (RAD23A antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.