

For Research Use Only

RAD23A Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86132-3-PBS



Basic Information

Catalog Number: 86132-3-PBS	GenBank Accession Number: BC014026	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 5886	CloneNo.: 250716F8
Source: Rabbit	UNIPROT ID: P54725	
Isotype: IgG	Full Name: RAD23 homolog A (S. cerevisiae)	
Immunogen Catalog Number: AG0716	Calculated MW: 40 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

86132-3-PBS targets RAD23A as part of a matched antibody pair:

MP02421-1: 86132-3-PBS capture and 86132-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



The mean RAD23A concentration was determined to be 190.9 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02421-1, RAD23A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86132-3-PBS. Detection antibody: 86132-2-PBS. Standard: Ag0716. Range: 0.312-40 ng/mL.



Bi-layer interferometry (BLI) kinetic assay of 86132-2-PBS against Human RAD23A was performed. The affinity constant is 0.113 nM.