

For Research Use Only

BTNL2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86137-3-PBS



Basic Information

Catalog Number: 86137-3-PBS	GenBank Accession Number: BC127642	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 56244	CloneNo.: 250562C10
Source: Rabbit	UNIPROT ID: Q9UIRO	
Isotype: IgG	Full Name: butyrophilin-like 2 (MHC class II associated)	
Immunogen Catalog Number: AG18803	Calculated MW: 455 aa, 50 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

86137-3-PBS targets BTNL2 as part of a matched antibody pair:

MP02518-1: 86137-2-PBS capture and 86137-3-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02518-1, Human BTNL2 Recombinant Matched Antibody Pair - PBS only. 86137-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag18803. 86137-3-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL



The mean BTNL2 concentration was determined to be 39.10 ng/mL in HeLa cell extract based on a 1.1 mg/mL extract load.