

For Research Use Only

REL Recombinant antibody, PBS Only (Capture)

Catalog Number: 86140-5-PBS



Basic Information

Catalog Number: 86140-5-PBS	GenBank Accession Number: NM_001291746	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 5966	CloneNo.: 250180B2
Source: Rabbit	UNIPROT ID: Q04864	
Isotype: IgG	Full Name: v-rel reticuloendotheliosis viral oncogene homolog (avian)	
Immunogen Catalog Number: AG26803	Calculated MW: 69 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

86140-5-PBS targets REL as part of a matched antibody pair:

MP02354-2: 86140-5-PBS capture and 86140-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean REL concentration was determined to be 43.7 ng/mL in Daudi cell extract based on a 1.1 mg/mL extract load, 3.9 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02354-2, REL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86140-5-PBS. Detection antibody: 86140-4-PBS. Standard: Ag26803. Range: 0.625-80 ng/mL.