

For Research Use Only

AKR1A1 Recombinant antibody, PBS Only

Catalog Number: 86147-1-PBS



Basic Information

Catalog Number:

86147-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7037

GenBank Accession Number:

BC000670

GeneID (NCBI):

10327

UNIPROT ID:

P14550

Full Name:

aldo-keto reductase family 1, member A1 (aldehyde reductase)

Calculated MW:

37 kDa

Observed MW:

36 kDa

Purification Method:

Protein A purification

CloneNo.:

250713E8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86147-1-RR (AKR1A1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86147-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 86147-1-RR (AKR1A1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86147-1-PBS in a different storage buffer formulation.