

For Research Use Only

GABARAP Recombinant monoclonal antibody

Catalog Number: 86165-3-RR



Basic Information

Catalog Number: 86165-3-RR	GenBank Accession Number: BC106748	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 11337	CloneNo.: 250354C7
Source: Rabbit	UNIPROT ID: O95166	Recommended Dilutions: WB: 1:5000-1:50000 IF/ICC: 1:250-1:1000
Isotype: IgG	Full Name: GABA(A) receptor-associated protein	
	Calculated MW: 14 kDa	
	Observed MW: 14-18 kDa	

Applications

Tested Applications: WB, IF/ICC, ELISA	Positive Controls: WB : HEK-293 cells, HeLa cells, SH-SY5Y cells, RAW 264.7 cells, NIH/3T3 cells, mouse brain tissue, rat brain tissue
Species Specificity: human, mouse, rat	IF/ICC : Chloroquine treated HepG2 cells,

Background Information

GABARAP (Gamma-aminobutyric acid type A receptor-associated protein), was initially identified as GABAA receptor binding protein and has shown to be implicated in intracellular trafficking of GABAA receptors. Through interacting with a broad range of proteins, GABARAP has been demonstrated as a multifunctional protein implicated in intracellular vesicle trafficking, apoptosis, and autophagy. The protein level of lipidated form of GABARAP (GABARAP-II) is upregulated during the autophagic process, which makes it as an autophagic marker. This antibody detected endogenous GABARAP protein in lysates from various tissues. It may cross-react with GABARAPL1 due to the high identity of amino acid sequence.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86165-3-RR (GABARAP antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Ethanol) fixed Chloroquine treated HepG2 cells using GABARAP antibody (86165-3-RR, Clone: 250354C7) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).