

For Research Use Only

Syntaxin 4 Recombinant antibody

Catalog Number: 86181-1-RR



Basic Information

Catalog Number:

86181-1-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6973

GenBank Accession Number:

BC002436

GeneID (NCBI):

6810

UNIPROT ID:

Q12846

Full Name:

syntaxin 4

Calculated MW:

34 kDa

Observed MW:

35-37 kDa

Purification Method:

Protein A purification

CloneNo.:

250825G3

Recommended Dilutions:

WB: 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : rat kidney tissue, mouse kidney tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 86181-1-RR (Syntaxin 4 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



rat kidney tissue were subjected to SDS PAGE followed by western blot with 86181-1-RR (Syntaxin 4 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.