

For Research Use Only

Frataxin Recombinant antibody

Catalog Number: 86195-4-RR



Basic Information

Catalog Number:

86195-4-RR

Size:

100ul, Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4622

GenBank Accession Number:

NM_000144.5

GeneID (NCBI):

2395

UNIPROT ID:

Q16595-1

Full Name:

frataxin

Calculated MW:

23kDa

Observed MW:

14-23 kDa

Purification Method:

Protein A purification

CloneNo.:

250761C8

Recommended Dilutions:

WB: 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : SW480 cells,

Background Information

Frataxin (FXN), also named as FRDA, X25, m81-FXN, d-FXN, m78-FXN and i-FXN, belongs to the frataxin family. It promotes the biosynthesis of heme and assembly and repair of iron-sulfur clusters by delivering Fe²⁺ to proteins involved in these pathways. FXN may play a role in the protection against iron-catalyzed oxidative stress through its ability to catalyze the oxidation of Fe²⁺ to Fe³⁺; the oligomeric form but not the monomeric form has in vitro ferroxidase activity. FXN is cleaved to be 4 chains. FXN is expressed in the cytoplasm as a 210 amino acid (AA) precursor protein (pFXN; 23 KDa) that is translocated into mitochondria where it is processed by two consecutive steps into iFXN (FXN 42-210; 19 KDa) and finally mFXN (81-210; 14.2 KDa), which is functional. (PMID: 26704351, PMID: 26671574)

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



SW480 cells were subjected to SDS PAGE followed by western blot with 86195-4-RR (Frataxin antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.