

For Research Use Only

alpha 2-Macroglobulin-like 1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86234-3-PBS



Basic Information

Catalog Number: 86234-3-PBS	GenBank Accession Number: GeneID (NCBI): 144568	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	UNIPROT ID: A8K2UO	CloneNo.: 250995C6
Source: Rabbit	Full Name: alpha-2-macroglobulin-like 1	
Isotype: IgG	Calculated MW: 161 kDa	
Immunogen Catalog Number: EG4523		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

86234-3-PBS targets alpha 2-Macroglobulin-like 1 as part of a matched antibody pair:

MP02375-1: 86234-3-PBS capture and 86234-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean A2ML1 concentration was determined to be 22.3 ng/mL in human urine.



Cytometric bead array standard curve of MP02375-1, A2ML1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86234-3-PBS. Detection antibody: 86234-2-PBS. Standard: Eg4523. Range: 1.562-200 ng/mL