

For Research Use Only

Draxin Recombinant antibody, PBS Only (Capture)

Catalog Number: 86244-2-PBS



Basic Information

Catalog Number: 86244-2-PBS	GenBank Accession Number: NM_198545.3	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 374946	CloneNo.: 250677C11
Source: Rabbit	UNIPROT ID: Q8NBI3	
Isotype: IgG	Full Name: chromosome 1 open reading frame 187	
Immunogen Catalog Number: EG4203	Calculated MW: 39 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

86244-2-PBS targets Draxin as part of a matched antibody pair:

MP02569-1: 86244-2-PBS capture and 86244-3-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02569-1, Human Draxin Recombinant Matched Antibody Pair - PBS only. 86244-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg4203. 86244-3-PBS was HRP conjugated as the detection antibody. Range: 0.781-50 ng/mL.



The mean Draxin concentration was determined to be 30.83 ng/mL in SH-SY5Y cell extract based on a 1.40 mg/mL extract load.