

For Research Use Only

CNPY3 Recombinant antibody

Catalog Number: 86256-2-RR



Basic Information

Catalog Number:

86256-2-RR

Size:

100ul, Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7382

GenBank Accession Number:

BC004423

GeneID (NCBI):

10695

UNIPROT ID:

Q9BT09

Full Name:

canopy 3 homolog (zebrafish)

Calculated MW:

31 kDa

Observed MW:

30-40 kDa

Purification Method:

Protein A purification

CloneNo.:

250998A1

Recommended Dilutions:

WB: 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HCT 116 cells, HEK-293 cells, Jurkat cells, THP-1 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86256-2-RR (CNPY3 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.