

For Research Use Only

SNRPA Recombinant antibody

Catalog Number: 86259-3-RR



Basic Information

Catalog Number: 86259-3-RR	GenBank Accession Number: BC000405	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 6626	CloneNo.: 251055G1
Source: Rabbit	UNIPROT ID: P09012	Recommended Dilutions: WB: 1:5000-1:50000 IF/ICC: 1:500-1:2000
Isotype: IgG	Full Name: small nuclear ribonucleoprotein polypeptide A	
Immunogen Catalog Number: AG0283	Calculated MW: 31 kDa	
	Observed MW: 31 kDa	

Applications

Tested Applications: WB, IF/ICC, ELISA	Positive Controls:
Species Specificity: human, mouse, rat	WB : MCF-7 cells, mouse brain tissue, HepG2 cells, HeLa cells, A431 cells, HEK-293 cells, rat thymus tissue IF/ICC : HepG2 cells,

Background Information

Heterogeneous nuclear ribonucleoproteins (hnRNPs) constitute a set of polypeptides that contribute to mRNA transcription, premRNA processing as well as mature mRNA transport to the cytoplasm and translation. They also bind heterogeneous nuclear RNA (hnRNA), which are the transcripts produced by RNA polymerase II. hnRNP A1 functions in the packaging of nascent RNA polymerase II transcripts and participates in a variety of nuclear and cytoplasmic processes that modulate gene expression. The RNA binding characteristics of hnRNP A1 suggest that it can modulate the expression of specific genes.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 86259-3-RR (SNRPA antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 86259-3-RR (SNRPA antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using SNRPA antibody (86259-3-RR, Clone: 251055G1) at dilution of 1:1000 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).



Biolayer interferometry (BLI) kinetic assays of 86259-3-RR against Human SNRPA were performed. The affinity constant is below 1 pM.