

For Research Use Only

Armetl1 Recombinant antibody

Catalog Number: 86264-1-RR



Basic Information

Catalog Number:

86264-1-RR

Size:

100ul , Concentration: 1000 µg/ml by 227526

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2845

GenBank Accession Number:

NM_177647.4

GeneID (NCBI):

227526

UNIPROT ID:

Q8CC36

Full Name:

arginine-rich, mutated in early stage tumors-like 1

Calculated MW:

21 kDa

Observed MW:

16-21 kDa

Purification Method:

Protein A purification

CloneNo.:

250599B3

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse

Positive Controls:

WB : mouse heart tissue, mouse testis tissue mouse skeletal muscle tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86264-1-RR (Armetl1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.