

For Research Use Only

TCF21 Recombinant antibody

Catalog Number: 86294-2-RR



Basic Information

Catalog Number:

86294-2-RR

Size:

100ul, Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG32510

GenBank Accession Number:

BC025697

GeneID (NCBI):

6943

UNIPROT ID:

O43680

Full Name:

transcription factor 21

Calculated MW:

179 aa, 20 kDa

Observed MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

250960C8

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : MG-63 cells,

Background Information

TCF21 (Transcription factor 21), also known as POD1. It is predicted to be located in the nucleus. TCF21 is a prognostic marker in Breast invasive carcinoma, Kidney renal clear cell carcinoma, Liver hepatocellular carcinoma, Lung squamous cell carcinoma. This protein is involved in the epithelial-mesenchymal interaction in the morphogenesis of kidney and lung, including epithelial differentiation and branching morphogenesis. It may play a role in the specialization or differentiation of one or more epicardial cell types. The molecular weight of TCF21 is 20 kDa, it also has phosphorylation modification.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86294-2-RR (TCF21 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 86294-2-RR against Human TCF21 were performed. The affinity constant is 0.797 nM.