

For Research Use Only

CRKRS Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86296-1-PBS



Basic Information

Catalog Number: 86296-1-PBS	GenBank Accession Number: BC150265	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 51755	CloneNo.: 250880C11
Source: Rabbit	UNIPROT ID: Q9NYV4	
Isotype: IgG	Full Name: Cdc2-related kinase, arginine/serine-rich	
Immunogen Catalog Number: AG25312	Calculated MW: 1490 aa, 164 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

86296-1-PBS targets CRKRS as part of a matched antibody pair.

MP02374-1: 86296-2-PBS capture and 86296-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CDK12/CRKRS concentration was determined to be 20.9 ng/mL in K-562 cell extract based on a 3.0 mg/mL extract load.



Cytometric bead array standard curve of MP02374-1, CDK12/CRKRS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86296-2-PBS. Detection antibody: 86296-1-PBS. Standard: Ag25312. Range: 0.781-100 ng/mL