

For Research Use Only

RPL17 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86348-1-PBS



Basic Information

Catalog Number: 86348-1-PBS	GenBank Accession Number: BC066323	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 6139	CloneNo.: 250901E3
Source: Rabbit	UNIPROT ID: P18621	
Isotype: IgG	Full Name: ribosomal protein L17	
Immunogen Catalog Number: AG5270	Calculated MW: 21 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

86348-1-PBS targets RPL17 as part of a matched antibody pair:

MP02355-1: 86348-2-PBS capture and 86348-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean RPL17 concentration was determined to be 17.6 ng/mL in HeLa cell extract based on a 1.5 mg/mL extract load.



Cytometric bead array standard curve of MP02355-1, RPL17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86348-2-PBS. Detection antibody: 86348-1-PBS. Standard: Ag5270. Range: 0.469-60 ng/mL.