

For Research Use Only

Cytokeratin 71 Recombinant antibody

Catalog Number: 86352-1-RR



Basic Information

Catalog Number:

86352-1-RR

Size:

100ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC103917

GeneID (NCBI):

112802

UNIPROT ID:

Q3SY84

Full Name:

keratin 71

Purification Method:

Protein A purification

CloneNo.:

250918C5

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : rat skin tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat skin tissue were subjected to SDS PAGE followed by western blot with 86352-1-RR (KRT71 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.