

For Research Use Only

MSRA Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86358-4-PBS



Basic Information

Catalog Number: 86358-4-PBS	GenBank Accession Number: BC054033	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 4482	CloneNo.: 251025C2
Source: Rabbit	UNIPROT ID: Q9UJ68	
Isotype: IgG	Full Name: methionine sulfoxide reductase A	
Immunogen Catalog Number: AG6053	Calculated MW: 26 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

86358-4-PBS targets MSRA as part of a matched antibody pair:

MP02504-1: 86358-4-PBS capture and 86358-3-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MSRA concentration was determined to be 1.82 ng/mL in HT-1080 cell extract based on a 1.4 mg/mL extract load.



Sandwich ELISA standard curve of MP02504-1, Human MSRA Recombinant Matched Antibody Pair - PBS only. 86358-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag6053. 86358-3-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL