

For Research Use Only

OXR1 Recombinant antibody, PBS Only (Detector)

Catalog Number:86359-1-PBS



Basic Information

Catalog Number: 86359-1-PBS	GenBank Accession Number: BC032710	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 55074	CloneNo.: 251027G6
Source: Rabbit	UNIPROT ID: Q8N573	
Isotype: IgG	Full Name: oxidation resistance 1	
Immunogen Catalog Number: AG4439	Calculated MW: 758 aa, 85 kDa	

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

86359-1-PBS targets OXR1 as part of a matched antibody pair:

MP02376-1: 86359-2-PBS capture and 86359-1-PBS detection (validated in Cytometric bead array, Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Sandwich ELISA standard curve of MP02376-1, Human OXR1 Recombinant Matched Antibody Pair - PBS only. 86359-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag4439. 86359-1-PBS was HRP conjugated as the detection antibody. Range: 78.1-2500 pg/mL



The mean OXR1 concentration was determined to be 21.4 ng/mL in K-562 cell extract based on a 3.0 mg/mL extract load.



The mean OXR1 concentration was determined to be 9,942.9 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load and 24,819.9 ng/mL in K-562 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02376-1, OXR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86359-2-PBS. Detection antibody: 86359-1-PBS. Standard: Ag4439. Range: 0.781-50 ng/mL