

For Research Use Only

Zyxin Recombinant antibody, PBS Only

Catalog Number: 86364-3-PBS



Basic Information

Catalog Number:

86364-3-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0419

GenBank Accession Number:

BC008743

GeneID (NCBI):

7791

UNIPROT ID:

Q15942

Full Name:

zyxin

Calculated MW:

80 kDa

Observed MW:

78 kDa

Purification Method:

Protein A purification

CloneNo.:

250924G5

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86364-3-RR (ZYX antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86364-3-PBS in a different storage buffer formulation.