

For Research Use Only

Zyxin Recombinant antibody

Catalog Number: 86364-3-RR



Basic Information

Catalog Number:

86364-3-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0419

GenBank Accession Number:

BC008743

GeneID (NCBI):

7791

UNIPROT ID:

Q15942

Full Name:

zyxin

Calculated MW:

80 kDa

Observed MW:

78 kDa

Purification Method:

Protein A purification

CloneNo.:

250924G5

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, HEK-293 cells, A431 cells, MCF-7 cells, Daudi cells, NIH/3T3 cells, PC-12 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1%BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86364-3-RR (ZYX antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.