

For Research Use Only

Mouse IL-25/IL-17E Recombinant monoclonal antibody, PBS Only

Catalog Number: 86396-1-PBS



Basic Information

Catalog Number:

86396-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3316

GenBank Accession Number:

NM_080729.3

GeneID (NCBI):

140806

UNIPROT ID:

Q8VHH8

Full Name:

interleukin 25

Calculated MW:

19KD

Observed MW:

20 kda

Purification Method:

Protein A purification

CloneNo.:

250033G5

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86396-1-RR (IL-25/IL-17E antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86396-1-PBS in a different storage buffer formulation.



rat testis tissue were subjected to SDS PAGE followed by western blot with 86396-1-RR (IL-25/IL-17E antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86396-1-PBS in a different storage buffer formulation.