

For Research Use Only

UFSP2 Recombinant antibody, PBS Only

Catalog Number: 86436-1-PBS



Basic Information

Catalog Number:

86436-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10692

GenBank Accession Number:

BC010493

GeneID (NCBI):

55325

UNIPROT ID:

Q9NUQ7

Full Name:

UFM1-specific peptidase 2

Calculated MW:

469 aa, 53 kDa

Observed MW:

53 kDa

Purification Method:

Protein A purification

CloneNo.:

251182F9

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, pig

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86436-1-RR (UFSP2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86436-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86436-1-RR against Human UFSP2 were performed. The affinity constant is below 1 pM.