

For Research Use Only

# XCL1 Recombinant antibody, PBS Only

Catalog Number: 86437-1-PBS



## Basic Information

Catalog Number:

86437-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2933

GenBank Accession Number:

NM\_002995.2

GeneID (NCBI):

6375

UNIPROT ID:

P47992

Full Name:

chemokine (C motif) ligand 1

Calculated MW:

13 kDa

Observed MW:

11 kDa

Purification Method:

Protein A purification

CloneNo.:

250232D1

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

XCL1 also known as Lymphotactin, is a member of the chemokine superfamily. Chemokines function in inflammatory and immunological responses, inducing leukocyte migration and activation.

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



PMA, Ionomycin and Brefeldin A treated NK-92 cells were subjected to SDS PAGE followed by western blot with 86437-1-RR (XCL1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86437-1-PBS in a different storage buffer formulation.