

For Research Use Only

ARD1A Recombinant antibody, PBS Only

Catalog Number: 86450-3-PBS



Basic Information

Catalog Number:

86450-3-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6519

GenBank Accession Number:

BC000308

GeneID (NCBI):

8260

UNIPROT ID:

P41227

Full Name:

ARD1 homolog A, N-acetyltransferase (S. cerevisiae)

Calculated MW:

26 kDa

Observed MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

251237F2

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86450-3-RR (ARD1A antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86450-3-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using ARD1A antibody (86450-3-RR, Clone: 251237F2) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 86450-3-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86450-3-RR against Human ARD1A were performed. The affinity constant is 5.33 nM.