

For Research Use Only

# NKG2C/CD159c Recombinant antibody, PBS Only (Capture)

Catalog Number: 86471-2-PBS



## Basic Information

<b>Catalog Number:</b> 86471-2-PBS	<b>GenBank Accession Number:</b> NM_002260.4	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug, Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 3822	<b>CloneNo.:</b> 250835A10
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P26717	
<b>Isotype:</b> IgG	<b>Full Name:</b> killer cell lectin-like receptor subfamily C, member 2	
<b>Immunogen Catalog Number:</b> EG1953	<b>Calculated MW:</b> 26kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Product Information

86471-2-PBS targets NKG2C/CD159c as part of a matched antibody pair:

MP02451-1: 86471-2-PBS capture and 86471-1-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP02451-1, Human NKG2C/CD159c Recombinant Matched Antibody Pair - PBS only. 86471-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg1953. 86471-1-PBS was HRP conjugated as the detection antibody. Range: 0.781-50 ng/mL



The mean NKG2C/CD159c concentration was determined to be 128.02 ng/mL in K-562 cell extract based on a 3.0 mg/mL extract load, 137.66 ng/mL in HL-60 cell extract based on a 3.7 mg/mL extract load and 150.96 ng/mL in HeLa cell extract based on a 1.4 mg/mL extract load.