

For Research Use Only

TPOR/MPL Recombinant monoclonal antibody, PBS Only

Catalog Number:86495-3-PBS



Basic Information

Catalog Number:

86495-3-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC153092

GeneID (NCBI):

4352

UNIPROT ID:

P40238

Full Name:

myeloproliferative leukemia virus oncogene

Calculated MW:

635 aa, 71 kDa

Observed MW:

84 kDa, 74 kDa

Purification Method:

Protein A purification

CloneNo.:

251731G4

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEL cells were subjected to SDS PAGE followed by western blot with 86495-3-RR (MPL antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86495-3-PBS in a different storage buffer formulation.



human peripheral blood platelets were subjected to SDS PAGE followed by western blot with 86495-3-RR (MPL antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86495-3-PBS in a different storage buffer formulation.