

For Research Use Only

PBX1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86504-1-PBS



Basic Information

Catalog Number: 86504-1-PBS	GenBank Accession Number: BC101578	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 5087	CloneNo.: 251373E1
Source: Rabbit	UNIPROT ID: P40424	
Isotype: IgG	Full Name: pre-B-cell leukemia homeobox 1	
Immunogen Catalog Number: AG12817	Calculated MW: 430 aa, 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

86504-1-PBS targets PBX1 as part of a matched antibody pair.

MP02469-1: 86504-2-PBS capture and 86504-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PBX1 concentration was determined to be 14.0 ng/mL in Hela cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02469-1, PBX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86504-2-PBS. Detection antibody: 86504-1-PBS. Standard: Ag12817. Range: 0.625-80 ng/mL.