

For Research Use Only

# HBXIP Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86515-2-PBS



## Basic Information

<b>Catalog Number:</b> 86515-2-PBS	<b>GenBank Accession Number:</b> BC062619	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 10542	<b>CloneNo.:</b> 251284H1
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> A0A8Z5A536	
<b>Isotype:</b> IgG	<b>Full Name:</b> hepatitis B virus x interacting protein	
<b>Immunogen Catalog Number:</b> AG5894	<b>Calculated MW:</b> 10 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

86515-2-PBS targets HBXIP as part of a matched antibody pair:

MP02485-1: 86515-2-PBS capture and 86515-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean HBXIP concentration was determined to be 348.1 ng/mL in K-562 cell extract based on a 1.2 mg/mL extract load, 334.9 ng/mL in U2OS cell extract based on a 1.3 mg/mL extract load.



Cytometric bead array standard curve of MP02485-1, HBXIP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86515-2-PBS. Detection antibody: 86515-1-PBS. Standard: Ag5894. Range: 0.625-80 ng/mL.