

For Research Use Only

LOXL1 Recombinant antibody

Catalog Number: 86534-1-RR



Basic Information

Catalog Number:

86534-1-RR

Size:

100ul, Concentration: 1000 µg/ml by

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24293

GenBank Accession Number:

BC015090

GeneID (NCBI):

4016

UNIPROT ID:

Q08397

Full Name:

lysyl oxidase-like 1

Calculated MW:

63 kDa

Observed MW:

60-63 kDa

Purification Method:

Protein A purification

CloneNo.:

251345F4

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : MCF-7 cells, MG-63 cells, PC-3 cells, U-87 MG cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86534-1-RR (LOXL1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.