

For Research Use Only

LAIR2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86559-3-PBS



Basic Information

Catalog Number: 86559-3-PBS	GenBank Accession Number: BC069366	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 3904	CloneNo.: 250338B11
Source: Rabbit	UNIPROT ID: Q6I554	
Isotype: IgG	Full Name: leukocyte-associated immunoglobulin-like receptor 2	
Immunogen Catalog Number: EG4247	Calculated MW: 152 aa, 16 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Product Information

86559-3-PBS targets LAIR2 as part of a matched antibody pair:

MP02633-1: 86559-3-PBS capture and 86559-2-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

