

For Research Use Only

RCAN2 Recombinant antibody, PBS Only

Catalog Number: 86561-1-PBS



Basic Information

Catalog Number: 86561-1-PBS	GenBank Accession Number: BC038509	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10231	CloneNo.: 251381D10
Source: Rabbit	UNIPROT ID: Q14206	
Isotype: IgG	Full Name: regulator of calcineurin 2	
Immunogen Catalog Number: AG3597	Calculated MW: 225 aa, 25 kDa	
	Observed MW: 27 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat, pig

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86561-1-RR (RCAN2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86561-1-PBS in a different storage buffer formulation.