

For Research Use Only

Mouse IL-17C Recombinant antibody, PBS Only

Catalog Number: 86563-1-PBS



Basic Information

Catalog Number: 86563-1-PBS	GenBank Accession Number: NM_145834.3	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 234836	CloneNo.: 250601E2
Source: Rabbit	UNIPROT ID: Q8K4C5	
Isotype: IgG	Full Name: interleukin 17C	
Immunogen Catalog Number: EG2840	Calculated MW: 22 kDa	

Applications

Tested Applications:
WB, Indirect ELISA
Species Specificity:
mouse


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein were subjected to SDS PAGE followed by western blot with 86563-1-RR (IL17c antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86563-1-PBS in a different storage buffer formulation.