

For Research Use Only

IL-17C Recombinant antibody

Catalog Number: 86563-1-RR



Basic Information

Catalog Number:

86563-1-RR

Size:

100ul , Concentration: 1000 µg/ml by

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2840

GenBank Accession Number:

NM_145834.3

GeneID (NCBI):

234836

UNIPROT ID:

Q8K4C5

Full Name:

interleukin 17C

Calculated MW:

22 kDa

Purification Method:

Protein A purification

CloneNo.:

250601E2

Recommended Dilutions:

WB: 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse

Positive Controls:

WB : Recombinant protein,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein were subjected to SDS PAGE followed by western blot with 86563-1-RR (IL17c antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.