

For Research Use Only

Mouse Il17c Recombinant antibody, PBS Only (Detector)

Catalog Number:86563-2-PBS



Basic Information

Catalog Number: 86563-2-PBS	GenBank Accession Number: NM_145834.3	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 234836	CloneNo.: 250601F7
Source: Rabbit	UNIPROT ID: Q8K4C5	
Isotype: IgG	Full Name: interleukin 17C	
Immunogen Catalog Number: EG2840	Calculated MW: 22 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
mouse

Product Information

86563-2-PBS targets Il17c as part of a matched antibody pair:

MP02646-1: 86563-3-PBS capture and 86563-2-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean IL-17C concentration was determined to be 3.61 ng/mL in mouse intestine tissue extract based on a 4.8 mg/mL extract load.



Sandwich ELISA standard curve of MP02646-1, Mouse IL-17C Recombinant Matched Antibody Pair - PBS only. 86563-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg2840. 86563-2-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL