

For Research Use Only

ZDHHC8 Recombinant monoclonal antibody, PBS Only

Catalog Number: 86572-2-PBS



Basic Information

Catalog Number:

86572-2-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26028

GenBank Accession Number:

BC053544

GeneID (NCBI):

29801

UNIPROT ID:

Q9ULC8

Full Name:

zinc finger, DHHC-type containing 8

Calculated MW:

83 kDa

Observed MW:

70-85 kDa

Purification Method:

Protein A purification

CloneNo.:

251438A7

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86572-2-RR (ZDHHC8 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86572-2-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86572-2-RR against Human ZDHHC8 were performed. The affinity constant is 0.915 nM.