

For Research Use Only

# ORC6L Recombinant monoclonal antibody, PBS Only

Catalog Number: 86591-1-PBS



## Basic Information

Catalog Number:

86591-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11769

GenBank Accession Number:

BC063565

GeneID (NCBI):

23594

UNIPROT ID:

Q9Y5N6

Full Name:

origin recognition complex, subunit 6 like (yeast)

Calculated MW:

252 aa, 28 kDa

Observed MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

251479G7

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)


W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86591-1-RR (ORC6L antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86591-1-PBS in a different storage buffer formulation.



mouse testis tissue were subjected to SDS PAGE followed by western blot with 86591-1-RR (ORC6L antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86591-1-PBS in a different storage buffer formulation.