

For Research Use Only

ORC6L Recombinant monoclonal antibody

Catalog Number: 86591-1-RR



Basic Information

Catalog Number: 86591-1-RR	GenBank Accession Number: BC063565	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 23594	CloneNo.: 251479G7
Source: Rabbit	UNIPROT ID: Q9Y5N6	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: origin recognition complex, subunit 6 like (yeast)	
Immunogen Catalog Number: AG11769	Calculated MW: 252 aa, 28 kDa	
	Observed MW: 28 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human, mouse

Positive Controls:

WB : HeLa cells, mouse testis tissue, Jurkat cells, MCF-7 cells, A549 cells, MDA-MB-231 cells

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86591-1-RR (ORC6L antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



mouse testis tissue were subjected to SDS PAGE followed by western blot with 86591-1-RR (ORC6L antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.