

For Research Use Only

PRDX2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86599-3-PBS



Basic Information

Catalog Number: 86599-3-PBS	GenBank Accession Number: BC003022	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 7001	CloneNo.: 251494G7
Source: Rabbit	UNIPROT ID: P32119	
Isotype: IgG	Full Name: peroxiredoxin 2	
Immunogen Catalog Number: AG0835	Calculated MW: 22 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Product Information

86599-3-PBS targets PRDX2 as part of a matched antibody pair:

MP02542-2: 86599-3-PBS capture and 86599-1-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



The mean Peroxiredoxin 2 concentration was determined to be 75.74 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load.



Plasma of six individual healthy human donors was measured. The human Peroxiredoxin 2 concentration of detected samples was determined to be 6.52 ng/mL with a range of 5.21 - 7.71 ng/mL.



Sandwich ELISA standard curve of MP02542-2, Human Peroxiredoxin 2 Recombinant Matched Antibody Pair - PBS only. 86599-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag0835. 86599-1-PBS was HRP conjugated as the detection antibody. Range: 3.13-200 ng/mL