

For Research Use Only

# MTHFD2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86608-2-PBS



## Basic Information

<b>Catalog Number:</b> 86608-2-PBS	<b>GenBank Accession Number:</b> BC017054	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 10797	<b>CloneNo.:</b> 251431A10
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P13995	
<b>Isotype:</b> IgG	<b>Full Name:</b> methylenetetrahydrofolate dehydrogenase (NADP+ dependent) 2, methenyltetrahydrofolate cyclohydrolase	
<b>Immunogen Catalog Number:</b> AG2911	<b>Calculated MW:</b> 350 aa, 38 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Product Information

86608-2-PBS targets MTHFD2 as part of a matched antibody pair:

MP02548-1: 86608-1-PBS capture and 86608-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS only, pH7.3


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean MTHFD2 concentration was determined to be 27.8 ng/mL in HeLa cell extract based on a 1.1 mg/mL extract load, 50.9 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02548-1, MTHFD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86608-1-PBS. Detection antibody: 86608-2-PBS. Standard: Ag2911. Range: 0.625-80 ng/mL.