

For Research Use Only

LRPPRC Recombinant antibody, PBS Only (Capture)

Catalog Number: 86616-1-PBS



Basic Information

Catalog Number: 86616-1-PBS	GenBank Accession Number: BC050311	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10128	CloneNo.: 251577C1
Source: Rabbit	UNIPROT ID: P42704	
Isotype: IgG	Full Name: leucine-rich PPR-motif containing	
Immunogen Catalog Number: AG15452	Calculated MW: 1394 aa, 158 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

86616-1-PBS targets LRPPRC as part of a matched antibody pair:

MP02559-1: 86616-1-PBS capture and 86616-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean LRPPRC concentration was determined to be 67.2 ng/mL in HeLa cell extract based on a 1.1 mg/mL extract load, 168.4 ng/mL in HCT 116 cell extract based on a 1.3 mg/mL extract load.



Cytometric bead array standard curve of MP02559-1, LRPPRC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86616-1-PBS. Detection antibody: 86616-2-PBS. Standard: Ag15452. Range: 1.562-100 ng/mL.