

For Research Use Only

INTS12 Recombinant antibody

Catalog Number: 86624-1-RR



Basic Information

Catalog Number: 86624-1-RR	GenBank Accession Number: BC014442	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 57117	CloneNo.: 251507A10
Source: Rabbit	UNIPROT ID: Q96CB8	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: integrator complex subunit 12	
Immunogen Catalog Number: AG9403	Calculated MW: 462 aa, 49 kDa	
	Observed MW: 49 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HepG2 cells, L02 cells
Species Specificity: human	

Background Information

INTS12 is a subunit of the Integrator complex that is thought to act as a reader to link the complex to chromatin.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86624-1-RR (INTS12 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.