

For Research Use Only

ALIX; AIP1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86644-4-PBS



Basic Information

| | | |
|---|--|---|
| Catalog Number: 86644-4-PBS | GenBank Accession Number: BC020066 | Purification Method: Protein A purification |
| Size: 100ug , Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 10015 | CloneNo.: 251574F8 |
| Source: Rabbit | UNIPROT ID: Q8WUM4 | |
| Isotype: IgG | Full Name: programmed cell death 6 interacting protein | |
| Immunogen Catalog Number: AG3074 | Calculated MW: 868 aa, 96 kDa | |

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Product Information

86644-4-PBS targets ALIX; AIP1 as part of a matched antibody pair:

MP02727-1: 86644-4-PBS capture and 86644-3-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



The mean Alix concentration was determined to be 44.88 ng/mL in HeLa cell extract based on a 1.20 mg/mL extract load.



Sandwich ELISA standard curve of MP02727-1, Human Alix Recombinant Matched Antibody Pair - PBS only. 86644-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag3074. 86644-3-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL.



Bi-layer interferometry (BLI) kinetic assay of 86644-4-PBS against Human ALIX; AIP1 was performed. The affinity constant is 1.62 nM.