

For Research Use Only

SAMD8 Recombinant monoclonal antibody, PBS Only

Catalog Number:86670-1-PBS



Basic Information

Catalog Number:

86670-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25310

GenBank Accession Number:

BC146556

GeneID (NCBI):

142891

UNIPROT ID:

Q96LT4

Full Name:

sterile alpha motif domain containing 8

Observed MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

251533B11

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86670-1-RR (SAMD8 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86670-1-PBS in a different storage buffer formulation.