

For Research Use Only

SAMD8 Recombinant monoclonal antibody

Catalog Number:86670-1-RR



Basic Information

Catalog Number: 86670-1-RR	GenBank Accession Number: BC146556	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 142891	CloneNo.: 251533B11
Source: Rabbit	UNIPROT ID: Q96LT4	Recommended Dilutions: WB: 1:1000-1:4000
Isotype: IgG	Full Name: sterile alpha motif domain containing 8	
Immunogen Catalog Number: AG25310	Observed MW: 45 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls:
Species Specificity: human, mouse, rat	WB : mouse heart tissue, rat heart tissue, human heart tissue

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86670-1-RR (SAMD8 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.