

For Research Use Only

SOX9 Recombinant monoclonal antibody

Catalog Number: 86714-3-RR



Basic Information

Catalog Number: 86714-3-RR	GenBank Accession Number: NM_000346	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 6662	CloneNo.: 253046H2
Source: Rabbit	UNIPROT ID: P48436	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: SRY (sex determining region Y)-box 9	
Immunogen Catalog Number: AG26274	Calculated MW: 56 kDa	
	Observed MW: 56 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : BxPC-3 cells, SW480 cells, HCT 116 cells, HT-29 cells, NCCIT cells, PC-3 cells
Species Specificity: human	

Background Information

SOX9, also named as Transcription factor SOX-9, is a 509 amino acid protein. Sex-determining region Y (SRY)-box 9 protein (SOX9) is a member of the SOX family of transcription factors (TFs) which are developmental regulators that possess high mobility group (HMG) box DNA binding and transactivation domains. It participates in a variety of functions, such as lineage restriction and terminal differentiation, through precise temporal and spatial expression patterns that differ between particular cell types and tissues. SOX9 gene has been implicated in different types of cancer as an oncogene; however, it also may behave as a tumor suppressor. SOX9 can be ubiquitinated or SUMOylated to higher molecular weight (PMID: 24155239; 16307912; 16554309; 32070068)

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86714-3-RR (SOX9 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 86714-3-RR against Human SOX9 were performed. The affinity constant is 1.12 nM.