

For Research Use Only

# TRIM47 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:86728-3-PBS



## Basic Information

<b>Catalog Number:</b> 86728-3-PBS	<b>GenBank Accession Number:</b> BC017299	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 91107	<b>CloneNo.:</b> 251436F2
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q96LD4	
<b>Isotype:</b> IgG	<b>Full Name:</b> tripartite motif-containing 47	
<b>Immunogen Catalog Number:</b> AG22868	<b>Calculated MW:</b> 638 aa, 70 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Product Information

86728-3-PBS targets TRIM47 as part of a matched antibody pair:

MP02626-1: 86728-3-PBS capture and 86728-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean TRIM47 concentration was determined to be 6.8 ng/mL in U2OS cell extract based on a 1.2 mg/mL extract load, 9.5 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 1.5 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02626-1, TRIM47 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86728-3-PBS. Detection antibody: 86728-2-PBS. Standard: Ag22868. Range: 0.312-40 ng/mL.