

For Research Use Only

NCF1/p47 phox Recombinant monoclonal antibody

Catalog Number:86733-3-RR



Basic Information

Catalog Number: 86733-3-RR	GenBank Accession Number: BC002816	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 653361	CloneNo.: 252310G3
Source: Rabbit	UNIPROT ID: P14598	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: neutrophil cytosolic factor 1	
Immunogen Catalog Number: AG28089	Calculated MW: 45 kDa	
	Observed MW: 45 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human, mouse

Positive Controls:

WB : Ramos cells, Daudi cells, Raji cells, RAW 264.7 cells

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86733-3-RR (NCF1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.