

For Research Use Only

MINA Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86760-1-PBS



Basic Information

Catalog Number: 86760-1-PBS	GenBank Accession Number: BC014928	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 84864	CloneNo.: 251692A8
Source: Rabbit	UNIPROT ID: Q8IUF8	
Isotype: IgG	Full Name: MYC induced nuclear antigen	
Immunogen Catalog Number: AG2838	Calculated MW: 465 aa, 53 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

86760-1-PBS targets MINA as part of a matched antibody pair:

MP02648-1: 86760-2-PBS capture and 86760-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



The mean MINA concentration was determined to be 23.7 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 29.1 ng/mL in A431 cell extract based on a 1.2 mg/mL extract load, 11.3 ng/mL in A375 cell extract based on a 1.6 mg/mL extract load.



Cytometric bead array standard curve of MP02648-1, MINA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86760-2-PBS. Detection antibody: 86760-1-PBS. Standard: Ag2838. Range: 0.469-60 ng/mL.



Biolayer Interferometry (BLI) kinetic assay of 86760-1-PBS against Human MINA was performed. The affinity constant is 3.19 nM.